



Therm-O-Flex® Appliances

INSTRUCTIONS AND CARE

General Instructions

1. Before insertion, soften appliance in warm water until it regains its flexibility
2. It is not necessary to store appliance in water when not in use
3. Clean appliance with a medium-soft tooth brush and nonabrasive tooth paste

NOTE: new appliances will become clear in approximately 1 week

Contact Us

If you have any questions or concerns, please don't hesitate to call us at 206-923-3684.

REFLECTIONS DENTAL CENTRE
2600 SW Barton Street, #E-20
Phone: (206) 923 – 3684
Fax: (206) 923 – 3368